LAB CAL

 [\(de\)en](#)

## Allgemein

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## General Information

The standard solutions for electrolytic conductivity have an exact conductivity value.

They are preferably used for the calibration of measuring instruments for electrolytic conductivity with highest accuracy.

Additional stabilizers prevent algae, bacteria and mould growth.

## Application

The standard solutions are ready to use.

The contents must be mixed before use in order to mix in any condensate in the upper part of the bottle. At temperatures below freezing point, segregation may occur with salt solutions. In such cases the solution should be tempered to 20 ... 25 °C and mix thoroughly.

It is state-of-the-art technology and good laboratory practice that the required quantity of standard is taken exclusively by pouring and of course is not poured back! This is the only way to ensure perfect, safe and traceable calibration of the measuring system.

Always keep the bottles tightly closed. Crystallization on the bottle or closure thread indicates that the bottles have not been sufficiently reclosed. Such bottles can no longer be used!

Residual contents of less than 25 % must be used immediately or discard.

## Calibration Instructions

- Shake bottle well before use.
- Remove and clean conductivity measuring cell.
- Rinse measuring cell and measuring vessel first with distilled water and then with reference solution. No dilution of the calibration solution must occur.
- Fill the measuring vessel with reference solution.
- For exact calibration: Adjust the temperature of reference solution exactly with a precision thermometer ( $\pm 0,05^\circ\text{C}$ ). to 25,0 °C.
- Follow the calibration steps as described in the instructions of the conductivity instrument used.

## Temperature tables

In most cases, the electrolytic conductivity is determined at exactly 25.0 °C (reference temperature) (laboratory) or uncalculated to this reference temperature using a table or temperature coefficients. We provide the following temperature tables for this purpose. Please contact us for temperature values that are not listed but are required.

## Calculations

The calculation of the temperature coefficient [TK %/K] is medium-specific and itself temperature-dependent. It must be determined for each measuring temperature for accurate measurements.

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## OEM/PLM Production

You are a device manufacturer and need calibration solutions with different conductivity values than here ? We manufacture every possible standard solution for you. Please contact us for further information.

## Notes

This product information exclusively relates to the product described in this leaflet. In particular, this product information cannot be applied to similar reagents from other manufacturers.

Periodically check for updates of this product information on our website.

## Instruction for use

For professional use only.

To avoid errors, the use of qualified personnel is carried out. National guidelines for work safety and quality assurance must be followed.

## Charts

The value tables listed here refer exclusively to standard solutions of Bioanalytic GmbH with the indicated article numbers. They cannot be used for products of other manufacturers.

## Storage

Store the conductivity standards at room temperature +15 ... +25 °C.

Always store conductivity standards away from gassing, strongly acidic or strongly alkaline solutions (e.g. hydrochloric acid, ammonia). Otherwise there is a risk of diffusion and a change in the conductivity values.

## Stability

In the original sealed condition, the standard solutions can be stored under the above storage conditions until the expiry date printed on them. Keep bottles tightly closed at all times.

Residual contents of less than 25 % bottle contents shall be used up immediately or discard.

During sampling, the removed quantity is replaced by ambient air and ambient air ( $\text{CO}_2$ ) is introduced. For this reason, standard conductivity solutions have only a limited shelf life after opening. Therefore, note the first opening date on the bottle.

## Support/Information service

For methodological and technical support, please contact us by E-Mail at [support@bioanalytic.de](mailto:support@bioanalytic.de) (German, English).

Periodically check for updates of this product information on our website.

## Feedback

Information from users can be reported to [support@bioanalytic.de](mailto:support@bioanalytic.de) (German, English).

Suggestions for further developments will be considered.

## Waste Management

Please observe your national laws and regulations.

Used and expired solutions must be disposed of in accordance with your local regulations. Inside the EU, national regulations apply that are based on the current, amended version of Council Directive 67/548/EEG on the approximation of the laws, regulations and administrative provisions relating to the classification, packaging and labelling of dangerous substances. Decontaminated packaging can be disposed of as household waste or recycled, unless otherwise specified.

## Unused Remains

These are usually hazardous wastes that must be recycled or disposed of. After consultation we take back such residual materials in the original container.

## Literature & Footnotes

Legends for the graphic symbols and tags used follow relevant norms or are available on our internet pages.

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## Electrolytic conductivity / temperature table

<b>Ganzzahlige Werte • Integer values</b>																		
°C:	041504		*) µS/cm	041505		041506		*) mS/cm	41507		041511		041512		04515 mS/cm	041513 mS/cm	041514 mS/cm	041516 mS/cm
	041504	*) µS/cm		041505	*) µS/cm	041506	*) µS/cm		41507	mS/cm	041511	mS/cm	041512	mS/cm				
0.0																		
5.0			316.3		633.0			0.6330		3.188					66.105	136.45	356.80	
10.0	71.8			360.8		721.0		0.7210		3.622					74.235	151.75	392.87	
15.0	80.9			406.5		808.8		0.8088		4.072					82.65	167.50	429.03	
20.0	90.3			453.0		904.8		0.9048		4.536					91.32	183.60	464.52	
<b>25.0</b>	<b>100.0</b>	<b>200.0</b>	<b>500.0</b>	<b>1000.0</b>				<b>1.0000</b>	<b>2.000</b>	<b>5.000</b>	<b>10.00</b>	<b>20.00</b>	<b>50.00</b>	<b>100.00</b>	<b>200.00</b>	<b>500.00</b>		
30.0	110.0			550.5		1099.3		1.0993		5.501					109.00	216.75	535.48	
35.0	120.2			602.5		1192.5		1.1925		6.001					118.10	233.60	568.47	
40.0	130.4			654.5		1301.0		1.3010		6.508					127.30	250.70	603.13	
45.0										7.024					136.50	136.50	636.05	
50.0				754.6		1510.3		1.5103		7.555					145.95	285.05	667.25	
55.0										8.102					155.45	155.45		
60.0										8.681					164.90	164.90		
65.0										9.215					174.50	174.50		
70.0										9.896					184.05			
75.0										10.510					194.05			
80.0										11.330					203.95			
85.0										12.165					214.00			
90.0																		
95.0																		

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