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Counting Chamber Nageotte

Maße:

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Metrics:

Chamber depth = 0.5 mm (special depths possible!)
Total area = 10 mm × 10 mm = 100 mm²; volume = 50 µL.
Division in 40 (rectangular) bands, each 10 × 0.25 mm side = 2.5 mm².
The horizontal ⁽²⁾ meridian is marked by a double line (one above and one below the median line).


Use:

Cell counting with very low cell values, e.g.:
CSF (Liquor cerebrospinalis), counting residual Leukocytes (WBC) in erythrocyte (RBC) and thrombocyte PLT concentrates or in plasma (transfusion medicine). Counting of nematodes.
The bands are counted by meandering succession

Calculation:

40 stripe fields: 2.5 mm² × 40 = 100 mm².
Chamber factor: 100 mm² × 0.5 mm = 50 µl. 1/5 = 0.02
Dilution factor: see your dilution
Calculation factor: 0.02 × dilution factor = calculation factor

Cell counts × calculation factor = cells/µL Sample material

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Calculation Examples

Leukocytes (WBCs)

Counting from WBC (white blood cells) and PLT (platelets) concentrates or from plasma. Count the entire chamber area in a meandering succession.

Dilution 1:20.

Standard dilution with Leuko-TIC® 1:20 *1).

Prepare by adding 20 µL sample to 380 µL Leuko-TIC®.

Calculation of calculation factor:

Chamber factor: $100\text{ mm}^2 \times 0,5\text{ mm} = 50\text{ }\mu\text{l}$. $1/50 = 0,02$
Dilution factor: $= 20$ (for dilution 1:20)
Calculation factor: $0,02 \times 20 = 0,4$

Calculation example:

Cell count in entire chamber area.

= 215

Multiply this result with the calculating factor **0.4**.

$215 \times 0,4 = 86/\mu\text{L}$ (= $0,086 \times 10^9/\mu\text{L} = 0,086 \times 10^9/\text{L} = 0,086$ Giga/L)

Dilution 1:5

Special dilution with Leuko-TIC®. Dilution 1:5 *1).

Prepare by adding 95 µL sample to 380 µL Leuko-TIC®.

Calculation of calculation factor:

Chamber factor: $100\text{ mm}^2 \times 0,5\text{ mm} = 50\text{ }\mu\text{l}$. $1/50 = 0,02$
Dilution factor: $= 5$ (for dilution 1:5)
Calculation factor: $0,02 \times 5 = 0,1$

Calculation example:

Cell count in entire chamber area.

= 215

Multiply this result with the calculating factor **0.1**.

$215 \times 0,1 = 21,5/\mu\text{L}$ (= $0,0215 \times 10^9/\mu\text{L} = 0,0215 \times 10^9/\text{L} = 0,0215$ Giga/L)

*1) Leuko-TIC® are single tests for quick, simple, clean and precise counting of WBCs. Also available: Ery-TIC®, Thombo-TIC®, Reti-TIC® and Sperm-TIC®

*2) Attention! Counting chambers with vertical lines are also available! These are for parts of the Asian market (vertical reading order). Routine counting takes place in a meandering structure from top left to bottom right.

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© Copyright by Bioanalytic GmbH